




500
pcs/case


3-layers


provide protection



Headrest Covers



HIGH
Quality

2 ply paper + 1ply poly

- ① PE waterproof protective film.
- ② Absorbent paper.
- ③ Absorbent paper.



3-LAYERS(2 PLY PAPER + 1PLY POLY)

The 3-layers design can provide optimal liquid barrier.

DISPOSABLE HEADREST COVERS:

Our softest comfort headrest cover with a pleasing appearance effectively protects patients and equipment.

CONVENIENT AND EASY TO USE:

Patented fabricel material provides absorbency and strength with a poly liner that prevents leak-through.

Reorder#	AB-HC-W
Size	10"×13" - 250mmx330mm
Color	White
Quantity	500 pieces/case



AB MED INC.

Address: 6021 Yonge Street, Unit 319
Toronto, Ontario, Canada M2M 3W2
Email: custserv@abmed.ca
Phone: +1 905-288-8199

